

COLD XENON LIGHT SOURCE

FIBROLUX X100

OPTIMUM ILLUMINATION FOR THE MEDICAL FIELD

- High intensity
- Electronic control
- Excellent performance



OPTOMIC[®]
producing quality

The Fibrolux X100 unit provides specialists of all medical fields with a high intensity light and excellent brightness for enhanced endoscopic visualization, using either flexible or rigid endoscopes. Likewise, it is ideal for microscopy and colposcopy.

COLD XENON LIGHT SOURCE **FIBROLUX X100**



DIFFERENT TYPES OF LIGHT GUIDE CONNECTORS



USER-FRIENDLY CONTROL PANEL



LIGHT AND PORTABLE

- Light power: 100W
- Excellent brightness, with light intensity adjustment
- Easy and quick lamp replacement
- Silent cooling system with automatic speed control
- Strong unit, manufactured in aluminium
- Simple operation
- Green or amber filter
- Power supply: 100 – 240V
- Dimensions: 150 x 120 x 280 mm
- Weight: 4.5 kg



Specifications are subject to change without notice

