

ENT RIGID ENDOSCOPES OTO-SINUS-LARYNGOSCOPES

OUTSTANDING QUALITY IN ENDOSCOPIC VISION


- High resolution
- Clear images
- High resistance



OPTOMIC offers a complete range of rigid otoscopes, sinusscopes and laryngoscopes for the inspection of the corresponding area. They allow for the diagnosis and serve as a method for minimal invasive surgical intervention. All the ENT OPTOMIC endoscopes are of proven quality and offer an exceptional image brightness to guarantee the satisfaction of the specialist.

ENT RIGID ENDOSCOPES


OTO-SINUS-LARYNGOSCOPIES




ROD LENSES SYSTEM



STAINLESS STEEL BODY



ADAPTERS



RESILIENT SAPPHIRE LENSES



PEEK EYEPIECE

0°

30°

45°

70°


70°

90°

	Ø	Length	0°	30°	45°	70°	70°	90°
OTOSCOPE	1,9	100	●	●				
	2,7	75	●	●		●	●	
	4	50	●	●				
SINUSCOPE	2,7	140	●	●		●	●	
	2,7	195	●	●		●	●	
	3	150	●	●	●	●		
	4	185	●	●	●	●		
LARYNGOSCOPE	4	185					●	
	6	158					●	
	10	158					●	●

● Non-Autoclavables ● Autoclavables

- Rod lenses system and stainless steel sealed by laser welding.
- Adapters for connection to type *Storz/Olympus*, *Wolf* and *ACMI* fiber optic cables.
- Autoclavable endoscopes with lenses made of sapphire of high resistance to scratches.



Specifications are subject to change without notice

